

肝細胞用培地 / Medium for Hepatocytes

製品番号	製品名	容量	保存温度 (配送方法)
肝細胞培養培地 Media and Plates for Hepatocyte Culture			
CM7500	Hepatocyte Thaw Medium (HTM)	45 mL	2 ~8 (冷蔵)
CM7000	Cryopreserved hepatocytes recovery medium (CHRM)	50 mL	2 ~8 (冷蔵)
CM3000	Primary Hepatocyte Thawing and Plating Supplements (Serum-containing)	for 500 mL medium	デキサメタゾン・カクテルA : 2 ~8 (冷蔵) FBS : -5 ~-20 (ドライアイス)
CM4000	Primary Hepatocyte Maintenance Supplements(Serum-free)	for 500mL medium	2 ~8 (冷蔵)
CM9000	Cryopreserved hepatocytes plating medium (CHPM)	50 mL	4 ~8 (冷蔵)
CM1096	Collagen-coated plate, 96-well	1 plate	15 ~30
A4000046701	Geltlex Flex LDEV-Free Reduced Glowth Factor Basement Membrane Matrix	1 mL	-20 ~-80 (ドライアイス)
A4000046702	Geltlex Flex LDEV-Free Reduced Glowth Factor Basement Membrane Matrix	5 mL	
A4000046703	Geltlex Flex LDEV-Free Reduced Glowth Factor Basement Membrane Matrix	10 mL	
A2737501	HepExtend™ Supplement (50X)	10 mL	-5 ~-20 (ドライアイス)
A1217601	WILLIAMS MED E W/O PHENOL RED 500ML	500 mL	2 ~8 (冷蔵)
12551032	WILLIAMS MED E	500 mL	2 ~8 (冷蔵)

供給元：サーモフィッシャーサイエンティフィック (ライフテクノロジーズ)

冷蔵禁止